

## Stereo Microscopes, Fixed Magnification, Two Powers



*ST-F802*

- Affordable
- Durable
- Easy to use

A basic instrument in many labs, classrooms and industrial work spaces, these work horse stereo microscopes feature:

- Clear, 3D images, erect, bright and flat
- Parfocal magnification stays in focus when switching magnification
- Illuminated stands-5W fluorescent below; 15W halogen above, or both
- Ample stage for sample manipulation yet fits in a storage cabinet
- Solid, metal focusing mechanisms endure usage beyond 5 year warranty

A special attraction to the ST-800 Series is the side mounted magnification control knobs, similar to a higher priced zoom stereo. Their all metal frame is built to last and low voltage bulbs last for 100+ hours.

Total Magnification with 10X eye piece	Jenco model no.	Field Diameter (mm)	Working Distance
10X / 20X	ST-F801	20.0 / 10.0	110mm.
10X / 30X	ST-F802	20.0 / 6.6	110mm.
10X / 40X	ST-F803	20.0 / 5.0	110mm.
20X / 40X	ST-F804	20.0 / 5.0	110mm.

*WHD 11 x 11 x 11*